

GRAMS STANDARD CANNULAE TIPS

MERCEDES TYPE

Three windows spaced equally around the tip.
2 Window also available.



#15	6mm x 32cm	5/8" HdL
#16	3.7mm x 30cm	5/8" HdL
#17	1.5mm x 8cm	1/2" HdL
#18	2.4mm x 14cm	1/2" HdL
#19	1.8mm x 14cm	1/2" HdL
#20	4.6mm x 30cm	5/8" HdL
#21	5.2mm x 30cm	5/8" HdL
#22	3mm x 23cm	5/8" HdL

FLAT TYPE 1

One Window Straight



#23	2.4mm x 14cm	1/2" HdL
#24	3mm x 14cm	1/2" HdL
#25	3.7mm x 14cm	1/2" HdL

FLAT TYPE 2

Two Slots in Line



#26	3.7mm x 14cm	5/8" HdL
#27	1.5mm x 10cm	1/2" HdL



SHARK TYPE

One window on the outside of the bend.
Tip Up

#4	6mm x 19cm	5/8" HdL
#5	5mm x 19cm	5/8" HdL
#6	3.7mm x 17cm	5/8" HdL
#7	6mm x 33cm	5/8" HdL
#9	3.7mm x 21cm	5/8" HdL
#10	2.4mm x 15cm	1/2" HdL
#11	6mm x 25cm	5/8" HdL
#13	3.7mm x 26cm	5/8" HdL



CONCORD TYPE

One window on the inside of the bend.
Tip Down

#8	6mm x 36cm	5/8" HdL
#12	6mm x 28cm	5/8" HdL
#14	3.7mm x 26cm	5/8" HdL

Custom Lengths Available

CANNULA HANDLE/HUB OPTIONS

5/8" Handle

1/2" Handle

60cc Toomey
Syringe Hub

SuperTwist LuerLock
Syringe Hub

- 5/8" x 5"
- All Stainless Steel
- Use 1/4", 3/8" silicone or 3/8" disposable tubing
- Use Part Number "MH5/8"
- 1/2" x 5"
- All Stainless Steel
- Use 1/4", 3/8" silicone or 3/8" disposable tubing
- Use Part Number "MH1/2"
- A hub connector for the Toomey Syringe
- For Syringe Aspirating
- All Stainless Steel
- Use Part Number "TH60" for Special Ordering
- A Superior LuerLock adapter
- All Stainless Steel
- Continuous thread design, engages the entire mating surfaces on the syringe threaded connection for maximum holding power to resist lateral and stripping forces.
- See standard mounted items featuring this connector
- Use Part Number "LH" for Special Ordering



Grams Medical

2443 Norse Ave – Costa Mesa, CA 92627
(949) 548-7337 – FAX (949) 548-1056

sales@gramsmmedical.com – www.gramsmmedical.com

SPECIALIZED INSTRUMENTS



GRAZER - GRAMS TIGER TIP® CANNULA

A more efficient instrument for the removal of tissue, particularly for revisions. This design is based on the Mercedes tip. The bullet nose is slotted back to the opening, giving the cannula an aggressive nature.


- Standard Diameters: 2.4, 2.7, 3, 3.7, 4.6, 5.2 mm.
- Lengths up to 40cm.
- Specify type of handle or hub, diameter and length when ordering.



HAMAS - WILLIAMS "GT" CANNULA *(Gynecomastia Tip)*

A very aggressive, specially modified Tiger Tip for Gynecomastia procedures. All leading edges are sharp.


- Standard Diameters: 3, 3.7, 4.6 mm.
- Standard Lengths: 14 & 20 cm.
- Available only on special 5.8" handle.
- Comes with protective silicone autoclavable tubing.



LAMBROS MULTI-HOLE CANNULA

A special tuned *finesse* cannula that has a controlled rate of removal. Produces smooth contours on soft tissue.

- Standard Diameter 3.0 mm.
- Standard Lengths: 24 & 30 cm. (custom lengths to 40 cm.)
- Specify 1/2" or 5/8" handle and length when ordering.



LAMBROS DIRECTIONAL CANNULA

Working with the multi-hole concept, but keeping all holes on a 180 plane makes this tip directional for superficial needs.

- Standard Diameter 3.0 mm.
- Standard Lengths: 24 & 30 cm. (custom lengths to 40 cm.)
- Specify 1/2" or 5/8" handle and length when ordering.



LAMBROS BUTTERFLY CANNULA

The Butterfly addition to the Lambros Cannula works as a spacer to preserve the subdermal tissue.


- Ideal for arms.
- Standard Diameter 3.4 mm.
- Custom lengths to 40 cm.
- Special 5/8" handle only.

NYLON BRUSHES

To clean cannula & handles

- A: 1.5 mm.
- B: 1.8 mm.
- C: 2.4 mm.
- D: 3 to 4.6 mm.
- E: 5.2 to 10 mm.

Use part # NB-"letter"
when ordering



Grams Medical

2443 Norse Ave – Costa Mesa, CA 92627
(949) 548-7337 – FAX (949) 548-1056

sales@gramsmedical.com – www.gramsmedical.com